



ORIGINAL RESEARCH

Open Access

# Midkine inhibitors: application of a simple assay procedure to screening of inhibitory compounds

Takashi Matsui<sup>1</sup>, Keiko Ichihara-Tanaka<sup>2,3</sup>, Chen Lan<sup>4</sup>, Hisako Muramatsu<sup>3,4</sup>, Toshiharu Kondou<sup>5</sup>, Chizuru Hirose<sup>5</sup>, Sadatoshi Sakuma<sup>1</sup> and Takashi Muramatsu\*<sup>2</sup>

## Abstract

**Background:** Midkine is a heparin-binding cytokine and is involved in etiology of various diseases. Thus, midkine inhibitors are expected to be helpful in treatment of many diseases.

**Methods:** We developed a simple assay for midkine activity based on midkine-dependent migration of osteoblastic cells. Midkine inhibitors were searched as materials that inhibit this midkine activity. To develop peptides that inhibit midkine activity, we constructed models in which C-terminal half of midkine interacted with  $\alpha_4\beta_1$ -integrin. Low molecular weight compounds which are expected to bind to midkine with high affinity were searched by *in silico* screening with the aid of Presto-X2 program.

**Results:** Among peptides in putative binding sites of midkine and the integrin, a peptide derived from  $\beta_1$ -integrin and that derived from the first  $\beta$  sheet of the C-terminal half of midkine significantly inhibited midkine activity. Two low molecular weight compounds found by *in silico* screening exhibited no toxicity to target cells, but inhibited midkine activity. They are trifluoro compounds: one (PubChem 4603792) is 2-(2,6-dimethylpiperidin-1-yl)-4-thiophen-2-yl-6-(trifluoromethyl)pyrimidine, and the other has a related structure.

**Conclusions:** The assay procedure is helpful in screening midkine inhibitors. All reagents described here might become mother material to develop clinically effective midkine inhibitors.

## Background

Midkine is a heparin-binding cytokine of molecular weight 13 kDa [1-3]. It enhances growth, survival and migration of various target cells. Midkine has around 50% sequence identity with pleiotrophin, and the two factors exhibit overlapping roles in many cases [1,4]. Midkine is also involved in initiation or progression of many pathological status, such as tumor invasion [5] rheumatoid arthritis [6], experimental autoimmune encephalitis [7], adhesion after surgery [8], neointima formation of the blood vessel [9], hypertension [10], and renal injury after ischemia [11], exposure to chemotherapeutic reagent [12] and diabetes [13]. Antisense oligonucleotides or siRNAs to midkine exhibit therapeutic effects in animal experiments concerning tumor growth [5,14-16], ischemic renal failure [17], neointima formation [18,19], adhesion

after surgery [20], and antibody-induced arthritis [20]. Polyclonal antibodies to midkine inhibit growth of Wilms' tumor cells *in vitro* [21]. Peptides or low molecular weight compounds are also expected to inhibit midkine activities, and in many cases these inhibitors might be superior to antisense oligoDNA or siRNA because of the ease in the administration methods.

In the present investigation, at first we developed a convenient method to screen midkine inhibitors. Then, applying the assay method, we tried to develop midkine inhibitors, namely peptides and other low molecular weight compounds.

Midkine is composed of two domains held by disulfide bridges. The C-terminal half domain is principally responsible for midkine activity and its heparin binding capability [22]. Two heparin binding sites are present in the C-terminal half [23]. Especially, Arg 81 has been identified to be important in heparin binding and midkine-induced neurite outgrowth and migration of neurons

\* Correspondence: tmurama@dpc.agu.ac.jp

<sup>2</sup> Department of Health Science, Faculty of Psychological and Physical Science, Aichi Gakuin University, 12 Arai-ke, Iwasaki-cho, Nisshin, Aichi 470-0195, Japan  
Full list of author information is available at the end of the article

[24,25]. Thus, we searched for peptides and low molecular weight compounds, which bind to the C-terminal half.

A rationale to obtain midkine inhibitory peptide should be based on midkine-receptor interaction. Midkine receptor is a molecular complex containing proteoglycans such as receptor-like protein tyrosine phosphatase  $\zeta$  [25], low density lipoprotein receptor-related protein [26] and integrin  $\alpha_4\beta_1$  or  $\alpha_6\beta_1$  [27]. A peptide sequence derived from low density lipoprotein receptor-related protein has already been utilized to develop an inhibitor of midkine [28]. In this study we were interested in interaction of midkine with  $\alpha_4\beta_1$  integrin.

## Methods

### Materials

Human midkine was produced in yeast as described before [29]. Candidate inhibitory peptides were synthesized by Peptide Institute, Osaka, Japan, as trifluoroacetate salts. Low molecular weight compounds for screening of midkine inhibitors were purchased from ChemDiv (San Diego, CA).

### Assay of midkine activity

Midkine activity was determined by promotion of migration of UMR106 (ATTC No. 1661) cells [30]. Chemotaxicell (pore size 8 micron, Kurabo Industries Ltd., Japan) was coated with 20  $\mu\text{g}/\text{ml}$  of midkine at the lower surface of the filter for 1 h, and washed with Dulbecco's phosphate-buffered saline (PBS) twice. UMR106 cells ( $2 \times 10^5$  cells in 0.2 ml of Dulbecco's modified Eagle's medium with 0.3% bovine serum albumin) were added to the upper chamber. The lower chamber contained 0.5 ml of the same medium. Cells were cultured for 4 h. The inner Chemotaxicell was cleaned with PBS, and cells were fixed with 100% methanol at room temperature for 20 min. The inner cell was cleaned with cotton bar, and the migrated cells were stained with 1% crystal violet at room temperature for 30 min., and were washed with  $\text{H}_2\text{O}$ . The cells were extracted with 0.2 ml of 1% SDS with 1% Triton X-100 for 1 h upon agitation. Then extract of 0.15 ml was transferred to a 96 well plate and OD 590 nm was determined. Upon inhibition assay using low molecular weight compounds, they were dissolved in ethanol and were added to the lower chamber so that the final ethanol concentration was 2%. The control contained also 2% ethanol in this occasion. Inhibition assay using peptides was performed by adding them dissolved in  $\text{H}_2\text{O}$ . Assay was done in duplicate.

### Other methods

Binding activity to midkine was evaluated in the BIA-CORE 3000 System (GE Healthcare, UK Ltd). A CM5 sensor chip coated with carboxymethyl dextran was activated by injecting a solution (1:1) of 200 mM N-ethyl-N'

(3-dimethylaminopropyl) carbodiimide hydrochloride and 100 mM N-hydroxysuccinimide at a flow rate of 10  $\mu\text{l}/\text{min}$ . A 70  $\mu\text{l}$  aliquot of midkine solution (50  $\mu\text{g}/\text{ml}$  in 10 mM acetate buffer, pH 5.0) was injected for immobilization. The remaining activated N-hydroxysuccinimide ester groups were blocked by injecting 1 M ethanolamine hydrochloride/NaOH, pH 8.5, and washed with 10  $\mu\text{l}$  of 1 M NaCl. Peptides (5  $\mu\text{g}/\text{ml}$ ) in the running buffer (10 mM HEPES-NaOH, pH 7.4, containing 0.15 M NaCl and 0.0005% Tween 20) were injected into the flow cell, and the change in resonance units was recorded. Binding assays were performed at 25°C with a constant flow rate of 10  $\mu\text{l}/\text{min}$  in both association and dissociation phases.

Docking of  $\alpha_4\beta_1$ -integrin with midkine was performed by the Greenpepper program (Chuden CTI) and docking of low molecular weight compounds in ZINC library <http://zinc.docking.org/> with midkine was performed by the Presto-X2 program (National Institute of Advanced Industrial Science and Technology, Tsukuba, Japan).

## Results

### Design of midkine inhibitory peptides based on interaction of midkine with integrin

The three dimensional structure of  $\alpha_4\beta_1$  integrin was deduced by homology modeling employing the structure of  $\alpha_v\beta_3$ -integrin [31]. The three dimensional structure of midkine was determined by NMR [23]. The possible three dimensional structures were deduced for the complex of the C-terminal half domain of midkine and  $\alpha_4\beta_1$ -integrin. Among 50 higher scores one, in 29 cases interactions occurred in relatively exposed area, namely, either the interacting site in the  $\alpha$ -chain is within 450 amino acid residues from the N-terminal and or interacting site in the  $\beta$ -chain is within 250 amino acid residues from the N-terminal. From the amino acid residues expected to be in the contact sites in the 29 cases, and from the homology of  $\alpha_v\beta_3$ - and  $\alpha_4\beta_1$ -integrins, we searched stretch of amino acids that have multiple contact sites and expected to be in  $\alpha$ -helix or  $\beta$ -sheet structure. The majority of contact sites in the integrin were in expected random coil structures, but sequences in Table 1 were in the expected  $\alpha$ -helix or  $\beta$ -sheet structures.

The expected binding sites in the midkine C-terminal half to  $\alpha_4\beta_1$  integrin were mostly in the three anti-parallel  $\beta$ -sheets [23]. Thus the following peptides in the  $\beta$ -sheets were selected as candidate peptide inhibitors derived from midkine: YKFENWGACDGG (YK12), GTGT-KVRQGTL (GT11), and QETIRVTKPC (QE10).

By surface plasmon assay, all of the integrin-derived peptides were found to have significant binding activity to midkine (Table 1). Then, we tested whether these peptides inhibited midkine activity.

**Table 1: Binding of integrin-derived peptides to midkine**

Name	Sequence	Origin	Expected structure	Binding constant to midkine ( $\mu\text{M}$ )
KQ9	KQNQVKFGS	$\alpha_4$ 247-255	$\alpha$ -helix	0.20
AM7	AMETNLV	$\alpha_4$ 359-365	$\beta$ -sheet	0.95
MQ7	MQSTIRE	$\alpha_4$ 335-341	$\beta$ -sheet	0.23
TF8	TFSQRIEG	$\alpha_4$ 419-426	$\beta$ -sheet	0.14
KA11	KANAKSCGECI	$\beta_1$ 29-39	$\alpha$ -helix	0.27
IE15	IENPRGSKDIKKKN	$\beta_1$ 80-94	$\beta$ -sheet	0.30
PA9	PAKLRNPCT	$\beta_1$ 200-208	$\alpha$ -helix	0.45

As the result YKFENWGACDGG (YK12) and KANAKSCGECI (KA11) were found to inhibit the migration activity when added at the concentration of 1 mg/ml to the lower chamber (Table 2). KA11 is derived from N-terminal region of  $\beta_1$ -integrin (Table 1). Midkine-derived YK12 is in an externally located  $\beta$ -sheet. When the peptides YK12 and KA11 were added to culture of UMR106 cells at the concentration of 1 mg/ml, no toxicity was observed, although inhibition of cell-substratum adhesion was observed in a certain occasion.


#### Search of low molecular weight compounds with inhibitory activity

The ZINC library of drug-like compounds were surveyed *in silico* for possible binding activity to C-terminal half of midkine. We selected 100 compounds, which were avail-

able from ChemDiv, one of the suppliers of drug-like compounds, and were among 531 compounds with the highest binding score in 300,000 compounds so far surveyed. Among the 100 compounds purchased, 10 compounds were soluble in ethanol at 1 mg/ml, and were not toxic to L cells after culturing for 24 h at the concentration of 40  $\mu\text{g}/\text{ml}$ . These 10 compounds were surveyed for the activity to inhibit midkine activity. Consequently, two compounds, Compound I [2-(2,6-dimethylpiperidin-1-

**Table 2: Inhibitory activities of synthetic peptides and low molecular weight compounds to midkine-dependent migration.**

Materials	Concentration	Migration activity (%)
None		100
Peptides		
YK12	1 mg/ml	44.0 $\pm$ 10.5
GT11	1 mg/ml	156 $\pm$ 6
QE10	1 mg/ml	99.5 $\pm$ 10.6
KQ9	1 mg/ml	107 $\pm$ 15
AM7	1 mg/ml	99.5 $\pm$ 6.9
MQ7	1 mg/ml	107 $\pm$ 2
TF8	1 mg/ml	105 $\pm$ 24
KA11	1 mg/ml	43.7 $\pm$ 12.8
IE15	1 mg/ml	113 $\pm$ 10.1
PA9	1 mg/ml	69.4 $\pm$ 46.1
Low molecular weight compounds		
Compound I	0.146 mg/ml	43.8 $\pm$ 10.2
Compound II	0.10 mg/ml	53.6 $\pm$ 4.1



yl)-4-thiophen-2-yl-6-(trifluoromethyl)pyrimidine; PubChem:4603792;ZINC00213280] and Compound II [ZINC00184372; no PubChem number, no official name] (Fig. 1) were found to significantly inhibit the activity (Table 2). They showed no toxicity to the target cells at the concentration used for the assay. The two compounds share a common structure, CF<sub>3</sub> (Fig. 1). In the model of the interaction of Compound I with the C-terminal half of midkine, CF<sub>3</sub> faces Lys 65 and Phe 66, which are in the first  $\beta$ -sheet (Fig 2). In the space filling model, compound I, the area of which is shown by a white line, also covered Arg. 81 (marked with an arrow as 1), which forms a key heparin-binding site [23] (Fig. 2). Mutation of Arg 81 leads to loss of neurite-promoting activity of midkine [24]. Compound II is also predicted to bind to midkine in a manner very similar to Compound I, while coverage of Arg 81 was only partial. (data not shown).

## Discussion

A convenient method has been developed to survey materials that inhibit midkine activity. This method has been effective to screen peptides and low molecular weight compounds with the inhibitory activity. An aptamer, which has midkine inhibitory activity and inhibits development of experimental autoimmune encephalitis, was also found using the method [7].

Concerning peptides and low molecular weight compounds with midkine binding and midkine inhibitory activity, these are expected to bind mainly to the first  $\beta$ -sheet of the C-terminal half domain. Probably, this strand is exposed and is easier to be bound by the inhibitors, and through this binding, R81, which is a key amino acid in midkine activity and comprises a heparin binding site [23,24], is expected to be covered by the inhibitor. In

addition, a possibility is not excluded that the first  $\beta$ -sheet itself plays a key role in midkine activity, such as in binding to integrins. Indeed a peptide derived from the first  $\beta$ -sheet was found to inhibit midkine activity. It is noted that the first  $\beta$ -sheet has Trp 69, which forms  $\beta$ C1 bulge with Arg81, the key residue in heparin-binding, which is in the centrally located  $\beta$ -sheet [23].

The inhibitory activity of peptides and low molecular weight compounds described here is not potent, and further elaboration is required to obtain materials of potential clinical utility. However, we think that the present material might even become seed materials. As an example, we noted that a methyl group in Compound I, a midkine inhibitory material, is expected to be located near to Arg 81 (marked with an arrow as 1 in Fig. 2), a key amino acid in midkine activity. Thus, conversion of the methyl group to carboxyl group might enhance the affinity of the compound to midkine, and increases the inhibitory activity.

## Conclusions

A simple assay procedure was developed for screening of midkine inhibitors. Inhibitory materials, especially low molecular weight compounds found by the present study are expected to become mother compounds to develop clinically applicable midkine inhibitors.

## Abbreviations

PBS: Dulbecco's phosphate-buffered saline

## Competing interests

The authors declare that they have no competing interests.

## Authors' contributions


TMa designed and performed an assay to screen midkine inhibitors, KI searched for low molecular weight compounds suitable for screening from candidate compounds, CL determined binding activity of synthetic peptides to midkine, HM designed candidate peptides with midkine inhibitory activity, TK designed the system of homology modeling and docking and performed *in silico* screening with CH, CH performed homology modeling, docking and *in silico* screening, SS designed over-all strategy to develop midkine inhibitors with TMu, and played key roles in preparing the second draft of the manuscript, TMu Designed over-all strategy to develop midkine inhibitors with S. S., designed *in silico* screening of low molecular weight compounds with midkine binding activity, and prepared the first draft of the manuscript. All authors contributed to the preparation of the manuscript and approved the final version.

## Acknowledgements

The Presto-X2 program was kindly made available by Dr. K. Fukunishi and Prof. H. Nakamura, National Institute of Advanced Industrial Science and Technology. This work was supported in parts by grants from the Ministry of Education, Culture, Sports, Science and technology of Japan.

## Author Details

<sup>1</sup>Cell Signals Inc, 3-29-18 Youkoudai, Isogoku, Yokohama 235-0045 Japan, <sup>2</sup>Department of Health Science, Faculty of Psychological and Physical Science, Aichi Gakuin University, 12 Araike, Iwasaki-cho, Nisshin, Aichi 470-0195, Japan, <sup>3</sup>Department of Nutritional Health, Faculty of Psychological and Physical Science, Aichi Gakuin University, 12 Araike, Iwasaki-cho, Nisshin, Aichi 470-0195, Japan, <sup>4</sup>Department of Biochemistry, Nagoya University Graduate School of Medicine, 65 Tsurumai-cho, Showa-ku, Nagoya 466-8550, Japan and <sup>5</sup>Chuden CTI Ltd, Higashisakura 1st Bld., 1-3-10 Higashisakura, Higashi-ku, Nagoya 461-0005, Japan



**Figure 2** A model of binding of Compound I to midkine as deduced by Presto X-2 program. White circle shows Compound I, and 1, 2 and 3 represent Arg81, Lys65 and Phe66, respectively.

Received: 20 April 2010 Accepted: 21 June 2010  
Published: 21 June 2010

## References

1. Muramatsu T: Midkine and pleiotrophin: two related proteins involved in development, survival, inflammation and tumorigenesis. *J Biochem* 2002, **132**:359-371.
2. Kadomatsu K, Muramatsu T: Midkine and pleiotrophin in neural development and cancer. *Cancer Lett* 2004, **204**:127-143.
3. Muramatsu T: Midkine, a heparin-binding cytokine with multiple roles in development, repair and diseases. *Proc Japan Acad Ser B* 2010, **86**:410-425.
4. Muramatsu H, Zou P, Kurosawa N, Ichihara-Tanaka K, Maruyama K, Inoh K, Sakai T, Chen L, Sato M, Muramatsu T: Female infertility in mice deficient in midkine and pleiotrophin, which form a distinct family of growth factors. *Genes Cells* 2006, **11**:1405-1417.
5. Takei Y, Kadomatsu K, Matsuo S, Itoh H, Nakazawa K, Kubota S, Muramatsu T: Antisense oligodeoxynucleotide targeted to Midkine, a heparin-binding growth factor, suppresses tumorigenicity of mouse rectal carcinoma cells. *Cancer Res* 2001, **61**:8486-8491.
6. Maruyama K, Muramatsu H, Ishiguro N, Muramatsu T: Midkine, a heparin-binding growth factor, is fundamentally involved in the pathogenesis of rheumatoid arthritis. *Arthritis Rheum* 2004, **50**:1420-1429.
7. Wang J, Takeuchi H, Sonobe Y, Jin S, Mizuno T, Miyakawa S, Fujiwara M, Nakamura Y, Kato T, Muramatsu H, Muramatsu T, Suzumura A: Inhibition of midkine alleviates experimental autoimmune encephalomyelitis through the expansion of regulatory T cell population. *Proc Natl Acad Sci USA* 2008, **105**:3915-3920.
8. Inoh K, Muramatsu H, Ochiai K, Torii S, Muramatsu T: Midkine, a heparin-binding cytokine, plays key roles in intraperitoneal adhesions. *Biochem Biophys Res Commun* 2004, **317**:108-113.
9. Horiba M, Kadomatsu K, Nakamura E, Muramatsu H, Ikematsu S, Sakuma S, Hayashi K, Yuzawa Y, Matsuo S, Kuzuya M, Kaname T, Hirai M, Saito H: Neointima formation in a restenosis model is suppressed in midkine-deficient mice. *J Clin Invest* 2000, **105**:489-495.
10. Hobo A, Yuzawa Y, Kosugi T, Kato N, Sato W, Maruyama S, Ito Y, Kobori H, Ikematsu S, Nishiyama A, Matsuo S, Kadomatsu K: The growth factor midkine regulates the rennin-angiotensin system in mice. *J Clin Invest* 2009, **119**:1616-1625.
11. Sato W, Kadomatsu K, Yuzawa Y, Hotta N, Matsuo S, Muramatsu T: Midkine is involved in neutrophil infiltration into the tubulointerstitium in ischemic renal injury. *J Immunol* 2001, **167**:3463-3469.
12. Kawai H, Sato W, Yuzawa Y, Kosugi T, Matsuo S, Takei Y, Kadomatsu K, Muramatsu T: Lack of the growth factor midkine enhances survival against cisplatin-induced renal damage. *Am J Pathol* 2004, **165**:1603-1612.
13. Kosugi T, Yuzawa Y, Sato W, Kawai H, Matsuo S, Takei Y, Muramatsu T, Kadomatsu K: Growth factor midkine is involved in the pathogenesis of diabetic nephropathy. *Am J Pathol* 2006, **168**:9-19.
14. Takei Y, Kadomatsu K, Goto T, Muramatsu T: Combinational antitumor effect of siRNA against midkine and paclitaxel on growth of human prostate cancer xenografts. *Cancer* 2006, **107**:864-873.
15. Dai LC, Wang X, Yao X, Lu YL, Ping JL, He JF: Enhanced therapeutic effects of combined chemotherapeutic drugs and midkine antisense oligonucleotides for hepatocellular carcinoma. *World J Gastroenterol* 2007, **13**:1989-1994.
16. Dai LC, Wang X, Yao X, Min LS, Ping JL, He JF: Antisense oligonucleotides targeting midkine inhibit tumor growth in an in situ human hepatocellular carcinoma model. *Acta Pharmacol Sin* 2007, **28**:453-458.
17. Sato W, Takei Y, Yuzawa Y, Matsuo S, Kadomatsu K, Muramatsu T: Midkine antisense oligodeoxyribonucleotide inhibits renal damage induced by ischemic reperfusion. *Kidney Int* 2005, **67**:1330-1339.
18. Hayashi K, Banno H, Kadomatsu K, Takei Y, Komori K, Muramatsu T: Antisense oligodeoxyribonucleotide as to the growth factor midkine suppresses neointima formation induced by balloon injury. *Am J Physiol Heart Circ Physiol* 2005, **288**:H2203-2209.
19. Banno H, Takei Y, Muramatsu T, Komori K, Kadomatsu K: Controlled release of small interfering RNA targeting midkine attenuates intimal hyperplasia in vein grafts. *J Vasc Surg* 2006, **44**:633-641.
20. Yamamoto H, Muramatsu H, Nakanishi T, Natori Y, Sakuma S, Ishiguro N, Muramatsu T: Midkine as a molecular target: comparison of effects of chondroitin sulfate E and siRNA. *Biochem Biophys Res Commun* 2006, **351**:915-919.
21. Muramatsu H, Shirahama H, Yonezawa S, Maruta H, Muramatsu T: Midkine, a retinoic acid-inducible growth/differentiation factor: immunochemical evidence for the function and distribution. *Dev Biol* 1993, **159**:392-402.
22. Muramatsu H, Inui T, Kimura T, Sakakibara S, Song X, Murata H, Muramatsu T: Localization of heparin-binding, neurite outgrowth and antigenic regions in midkine molecule. *Biochem Biophys Res Commun* 1994, **203**:1131-1139.
23. Iwasaki W, Nagata K, Hatanaka H, Inui T, Kimura T, Muramatsu T, Tasumi M, Inagaki F: Solution structure of midkine, a new heparin-binding growth factor. *EMBO J* 1997, **16**:6936-6946.
24. Asai T, Watanabe K, Ichihara-Tanaka K, Kaneda N, Kojima S, Iguchi A, Inagaki F, Muramatsu T: Identification of heparin-binding sites in midkine and their role in neurite-promotion. *Biochem Biophys Res Commun* 1997, **236**:66-70.
25. Maeda N, Ichihara-Tanaka K, Kimura T, Kadomatsu K, Muramatsu T, Noda M: A receptor-like protein-tyrosine phosphatase PTP $\zeta$ /RPTP $\beta$  binds a heparin-binding growth factor midkine. Involvement of arginine 78 of midkine in the high affinity binding to PTP $\zeta$ . *J Biol Chem* 1999, **274**:12474-12479.
26. Muramatsu H, Zou K, Sakaguchi N, Ikematsu S, Sakuma S, Muramatsu T: LDL receptor-related protein as a component of the midkine receptor. *Biochem Biophys Res Commun* 2000, **270**:936-941.
27. Muramatsu H, Zou P, Suzuki H, Oda Y, Chen GY, Sakaguchi N, Sakuma S, Maeda N, Noda M, Takada Y, Muramatsu T:  $\alpha$ -4 $\beta$ 1 and  $\alpha$ -6 $\beta$ 1 integrins are functional receptors for midkine, a heparin-binding growth factor. *J Cell Sci* 2004, **117**:5405-5415.
28. Chen S, Bu G, Takei Y, Sakamoto K, Ikematsu S, Muramatsu T: Midkine and LDL-receptor-related protein 1 contribute to the anchorage-independent cell growth of cancer cells. *J Cell Sci* 2007, **120**:4009-4015.
29. Ikematsu S, Yano A, Aridome K, Kikuchi M, Kumai H, Nagano H, Okamoto K, Oda M, Sakuma S, Aikou T, Muramatsu H, Kadomatsu K, Muramatsu T: Serum midkine levels are increased in patients with various types of carcinomas. *Br J Cancer* 2000, **83**:701-706.
30. Qi M, Ikematsu S, Maeda N, Ichihara-Tanaka K, Sakuma S, Noda M, Muramatsu T, Kadomatsu K: Haptotactic migration induced by midkine. Involvement of protein-tyrosine phosphatase  $\zeta$ , mitogen-activated protein kinase, and phosphatidylinositol 3-kinase. *J Biol Chem* 2001, **276**:15868-15875.
31. Xiong JP, Stehle T, Diefenbach B, Zhang R, Dunker R, Scott DL, Joachimiak A, Goodman SL, Arnaout MA: Crystal structure of the extracellular segment of integrin  $\alpha$ V $\beta$ 3. *Science* 2001, **294**:339-345.

doi: 10.1186/1755-7682-3-12

Cite this article as: Matsui et al., Midkine inhibitors: application of a simple assay procedure to screening of inhibitory compounds *International Archives of Medicine* 2010, **3**:12

Submit your next manuscript to BioMed Central and take full advantage of:

- Convenient online submission
- Thorough peer review
- No space constraints or color figure charges
- Immediate publication on acceptance
- Inclusion in PubMed, CAS, Scopus and Google Scholar
- Research which is freely available for redistribution

Submit your manuscript at  
www.biomedcentral.com/submit

